

WALL MOUNTED HEALTHCARE BIN RANGE

CODE: W2080


WALL MOUNTED LARGE HEALTHCARE BIN



Specifications

- ◆ **Wall Mounted Bin**
- ◆ **Space to Clean under Bin**
- ◆ **Stainless Steel Pedal**
- ◆ **Anti-bacterial Powdercoat Finish**
- ◆ **Screening or Labelling**
- ◆ **Silent Closing Option**
- ◆ **Heavy Duty Locking System**
- ◆ **Medical Grade Non Porous Shock Cord**


FITTING INSTRUCTIONS



Stage 1

Locate Bracket on Wall at a Height of 718mm from Floor to Bottom of Bracket

Fix using Toggle Bolts or similar Bolts suitable use on Plasterboard



Stage 2

Special Fixing stand is supplied with All Wall Mounted Bins Place on floor under Bracket



Stage 3

Place Bin on top of Special Fixing Stand & move closer to the Wall



Stage 4

Ensure Fixing Stand is flush with Front of Pedal



Stage 5

Move Bin into position



Stage 6

Attached Bin to Bracket using Supplied Bolts using M4 Allen Key & remove fixing stand

ALL PRODUCTS DESIGNED & MANUFACTURED IN IRELAND

COPYRIGHT © 2004

Paragon
PRODUCTS

TALTECH ENGINEERING
HOLDINGS LTD.
Unit 4 / Unit 5
Jamestown Rd. Ind. Est.
Inchicore, Dublin 8.

Tel: 01 453 6600
Fax: 01 478 3069
Email: info@taltech.ie
www.taltech.ie
www.paragonproducts.ie

